### Hardness of water

# Water (H2O)

## Soft water

## Hard water

#### Soft water

- Water consists of low concentration of calcium and magnesium salts.
- It gives foam with soap
- Examples:
  - Tap water.
  - Drinking water.

#### Hard water

- Hard water is due to the presence of high concentration of calcium and magnesium salts that are dissolved in water.
- It doesn't form foam with soap.
- Examples:
  - Sea water.

### Hard water

Temporary

Hard water

Permanent Hard water

#### **Temporary Hardness**

- Temporary Hardness is due to the presence of Ca<sup>2+</sup>, Mg<sup>2+</sup> in the form of the bicarbonate ion HCO<sub>3</sub><sup>-</sup>, being present in the water.
- This type of hardness can be treated by boiling the water to expel the CO<sub>2</sub>, as indicated by the following equation:
- $Ca^{2+} + 2HCO_3^{-} \rightarrow CaCO_3 + H_2O + CO_2$
- $Mg^{2+} + 2HCO_3^- \rightarrow MgCO_3 + H_2O + CO_2$
- Bicarbonate hardness is classified as temporary hardness

#### **Permanent Hardness**

- Permanent hardness is due to the presence of the ions Ca<sup>2+</sup>, Mg<sup>+2</sup> in the form of Cl<sup>-</sup>and SO<sub>4</sub><sup>2-</sup>. This type of hardness cannot be eliminated by boiling.
- The water with this type of hardness is said to be *permanently hard*
- As it can't be treated easily ,so it's treated by chemical treatment such as: ion exchange resin.

#### **Problems of hard water**

- Originally, water hardness was defined as the measure of the capacity of the water to precipitate soap
- It forms scales in the boiler that may cause:
  - 1. Decreasing in heat exchange efficiency.
  - 2. Corrosion takes place.
  - 3. Explosion.

#### **Scales due to hard water**



#### **Scales due to hard water**



Boiler scale on water side

#### **Determination of Total hardness**

The ions involved in water hardness, i.e. Ca<sup>2+</sup>(aq) and Mg<sup>2+</sup>(aq), can be determined by titration with a chelating agent ethylenediaminetetraacetic acid (EDTA), usually in the form of disodium salt (H<sub>2</sub>Y<sub>2</sub>-). The titration reaction is:

During titration: 
$$\begin{array}{cccc} H_2Y^{2-}(aq) &+ & Ca^{2+}(aq) \\ H_2Y^{2-}(aq) &+ & Mg^{2+}(aq) \end{array} \xrightarrow{\rightarrow} & CaY^{2-}(aq) &+ & 2 H^+(aq) \\ H_2Y^{2-}(aq) &+ & Mg^{2+}(aq) \end{array} \xrightarrow{\rightarrow} & MgY^{2-}(aq) &+ & 2 H^+(aq) \\ \end{array}$$
At end point 
$$\begin{array}{cccc} H_2Y^{2-}(aq) &+ & MgIn^-(aq) \\ Wine-red \end{array} \xrightarrow{\rightarrow} & MgY^{2-}(aq) &+ & HIn^{2-}(aq) \\ Wine-red & & sky-blue \end{array}$$

#### **Determination of Total hardness**



### **Determination of Total hardness**

- When both Ca and Mg are both determined, this experiment is called total hardness.
- **Chelating agent:** EDTA which is capable to react with Ca and Mg ions that present in sample solution.
- Hardness solution(I), (II): (NH<sub>4</sub>OH+NH<sub>4</sub>Cl) and (KOH) used as buffer solution that keeps the pH = 10.
- Sample solution: Tap water.
- Indicators: ManVer which has the ability to select both (Ca, Mg) ions and the calVer has the ability to select Ca ion only.

### **Determination of total hardness**

#### **Procedures:**

- Get 10 ml of tap water as a sample solution in a conical flask.
- Add 1 ml of hardness solution (I) to the sample.
- Add 2 drops of indicator (manVer) to the sample solution.
- Fill the burette with standard solution of EDTA (0.01M).
- Titrate EDTA against the sample until the color of the indicator changes from red to blue.
- Repeat these steps 3 times.
- Calculate the average value of the three volumes.

#### **Determination of total hardness**

Reading	V <sub>1</sub> (ml)	<b>V</b> <sub>2</sub> (ml)	<b>V</b> <sub>3</sub> (ml)	V <sub>average</sub> (ml)
Volume	•••••	•••••	•••••	•••••

$$V_{average} = (V_1 + V_2 + V_3) / 3$$

Concentration of  $(Ca + Mg) = V_{av} * 100 = \dots mg/I$ 

#### **Determination of calcium hardness**

#### **Procedures:**

- Get 10 ml of tap water as a sample solution in a conical flask.
- Add 1 ml of KOH as hardness solution to the sample.
- Add 2 drops of indicator (calVer) to the sample solution.
- Fill the burette with standard solution of EDTA (0.01M).
- Titrate EDTA against the sample until the color of the indicator changes from red to blue.
- Repeat these steps 3 times.
- Calculate the average value of the three volumes.

#### **Determination of calcium hardness**

Reading	V <sub>1</sub> (ml)	<b>V<sub>2</sub> (</b> ml)	<b>V</b> <sub>3</sub> (ml)	V <sub>average</sub> (ml)
Volume	•••••	•••••	•••••	

$$V_{average} = (V_1 + V_2 + V_3) / 3$$

Concentration of (Ca) =  $V_{av}$  \*100 = ..... mg/l

